

ELSODENT

Clinical and more biocompatible products

Elsotemp

Eugenol-free, fluoridated temporary filling cement



Features

- Because ELSOTEMP contains a high quality plasticizer, no expansion will occur during the setting process.
- No pulp or gum irritations (PH=7).
- Fast-setting.
- Waterproof.
- Eugenol-free: does not hinder polymerization of resins and composites.
- Contains fluoride: cariostatic action.
- Soft before setting to allow a perfect application during the filling process.
- Hard after setting.
- Stick to the walls of the tooth: ELSOTEMP can remain in mouth 1 or 2 months, if needed.
- Easy to remove.

Consistency

- Grey: soft.
- White: normal.
- Pink: hard.

References & Presentations

TPS-B	Temporary filling cement, white, 3 x 28 g syringes.
TPS-R	Temporary filling cement, pink, 3 x 28 g syringes.
TPP-R	Temporary filling cement, white, 6 x 28 g jars.
TPP-B	Temporary filling cement, pink, 6 x 28 g jars.
TPP-G	Temporary filling cement, grey, 6 x 28 g jars.
TPP-R1	Temporary filling cement, pink, 1 x 28 g jar.
TPP-B1	Temporary filling cement, white, 1 x 28 g jar.
TPP-G1	Temporary filling cement, grey, 1 x 28 g jar.
TPP-R40	Temporary filling cement, pink, 1 x 40 g jar.
TPB-R40	Temporary filling cement, white, 1 x 40 g jar.

MD: class IIA. Certifying body: SNCH. Manufacturer: G-Pharma



Phone: +33 (0)1 30 37 75 75 - elsodent@elsodent.com www.elsodent.com
ELSODENT, 35 Avenue du Gros Chêne - BP 50072 - HERBLAY - 95612 CERGY PONTOISE CEDEX - FRANCE