

ELSODENT

Clinical and more biocompatible products

Hemosal

Astringent gel – 25% aluminum sulfate



Features

- Thixotropic gel: easy, precise and fast setting.
- Due to its balance viscosity, the gel remains in place during all the procedure
- No tissue necrosis or blackening of tissue.
- Versatile in uses.
- Blue colour for easy visibility and placement.
- No contra-indications.

Indications

- In restorative dentistry: conditioning of the injured gum during the preparation of juxta-gingival class V cavities, to avoid discoloration of the definitive composite restoration.
- Prosthetic preparation: to allow a non-traumatic and optimal retraction of the gum. Apply the gel in the sulcus and insert the retraction cord.
- In endodontic procedures: when a canal is bleeding.



References & Presentations

HSA-6	4 x 1.5g syringes + 10 black needle tips.
LTL	100 intra-oral black tips.

MD: class IIA. Certifying body: SNCH. Manufacturer: G-Pharma



Phone: +33 (0)1 30 37 75 75 - elsodent@elsodent.com www.elsodent.com
ELSODENT, 35 Avenue du Gros Chêne - BP 50072 - HERBLAY - 95612 CERGY PONTOISE CEDEX - FRANCE