

## CORE D FLOW

Matériau composite «dual» de reconstitution


Sans HEMA, TEGDMA

### Points forts

**Thixotropie étudiée**, pour l'injection dans les cavités difficiles d'accès et quand la fluidité n'est pas un problème pour la mise en place.

Scellement des tenons radiculaires, en conjonction avec un adhésif «dual».

Taux de conversion, 68%.



Excellente résistance à la flexion, 136 MPa.

**Dureté comparable à celle de la dentine**, aucune différence de sensation quand la fraise passe de la dentine au matériau.

### Indications

- Reconstitution composites (dents dévitalisées, etc.)
- Scelllements des tenons radiculaires (en conjonction avec HEALBOND MP)

### Références & Présentations

<b>CDF-5</b> (dentine)	2 seringues dual de 2,5ml + 10 embouts mélangeurs marrons + 10 embouts intra-oraux.
<b>EMP-100</b>	Embouts mélangeurs, bagues marron. quantité 100.
<b>IMP-100</b>	100 embouts intra-oraux médiums.


Dispositif médical de classe IIa certifié par Tüv Rheinland (0197).

## CORE D FLOW

Dual-curing core build-up material

HEMA and TEGDMA -free

### Features



### Indications

- Core built-up on non-vital teeth.
- As a composite on root canal around the posts.

### References & Presentations

**CDF-5** (dentin) 2 x 2,5ml syringes + 10 mixing tips + 5 intra-oral regular tips + 5 intra-oral x-fine tips.

MD : Class IIa certified by Tüv Rheinland (0197).